

Day 10: Design 3

Sunday, July 30, 2023 3:21 PM

Goals:

- Continue working on software design
 - Pseudocode for app and RTC init still need to be written
 - Revising DAQ pseudocode based on preliminary testing
- Begin developing some of the software for the DAQ and App
- Aid with other parts of the design, time permitting

RTC initialization:

Arduino Code responsible to set up the RTC, runs once.

External Libraries

Arduino
uRTCLib
UnixTime
WiFinINA

Globals:

ssid: char[]
 The WiFi SSID
pass: char[]
 The WiFi password
rtc: uRTCLib::uRTCLib
 uRTCLib library object
stamp: UnixTime::UnixTime
 Unix timestamp converter
status: int
 WiFi radio status

Arduino Functions:

setup: () => void

 Begin RTC module

 Set status to WiFinINA::WL_IDLE_STATUS

 While the status is not WiFinINA::WL_CONNECTED:
 Attempt to connect (Begin WiFi, status is updated)
 Wait for 10 seconds

 Set RTC epoch using Wifi.getTime()

loop: () => void

 EMPTY

App Code:

Framework: Ionic/React

Using "Tabs" template as a baseline

External Libraries:

- @ionic/react
- @capacitor-community/bluetooth-le
- @capacitor/core
- @capacitor/filesystem
- @capacitor/preferences
- @capacitor/local-notifications
- ionic-simple-lockscreens
- ionicons
- chart.js
- react-chartjs-2

Hooks:

Globals:

- HRV_SERVICE: const string
Bluetooth Low Energy HRV Metric Service UUID
- HRV_CHARACTERISTIC: const string
Bluetooth Low Energy Characteristic UUID (receive HRV records)
- ERROR_CHARACTERISTIC: const string
Bluetooth Low Energy Characteristic UUID (receive error codes from wearable).
- REQUEST_CHARACTERISTIC: const string
Bluetooth Low Energy Characteristic UUID (send data requests to wearable).

data_hook: async (callbacks: array) => void

Try:

- Initialize the bluetooth-le BleClient
- Let deviceId: string be the output of connect(callbacks)

If there is an error:

- Log the error in console

connect: async (callbacks) => string

Let connected:boolean be false

While not connected:

- If there is no saved id in preferences:
 - Search for a device with the HRV_SERVICE as its advertised service

- Save selected device id to preferences

Else:

- Use saved id from preferences

- Attempt to connect to the device wearable 5 times with a timeout of 5 seconds
(On Disconnect parameter is connect(callbacks) function).

If no connection:

- Clear preference for id

Let lastWeekDate: string be the UTC ISO8601 representation of 1 week ago

Write lastWeekDate to REQUEST_CHARACTERISTIC.

Start Notifications for HRV and Error characteristics

On Disconnect Parameters:

HRV_CHARACTERISTIC => callbacks[0]
ERROR_CHARACTERISTIC => callbacks[1]

fetchRecords: async (timePeriod:number) => array
Let rawData: array be an empty array

Let endingTimestamp: number be the current unix timestamp

Let currentDatetime: Date be a UTC Date object initialized to endingTimestamp - timePeriod

Let startTimestamp be currentDatetime's unix timestamp

While the unix time of currentDatetime is before endingTimestamp:
Try to open the corresponding HRV data file using currentDatetime (req...)

If the file doesn't exist
Increment currentDatetime by a day
Continue with the loop

Process data file:
Split each record (separated by a newline)
Split the numbers within each record (separated by a space).
Add each element from resulting array from split to rawData

Filter rawData by removing records with timestamps lower than startingTimestamp and greater than endingTimestamp

Return rawData

Pages:

Home Screen:

HomeScreen:React.FC: ({userState}) => JSX

Return JSX object implementing the design of ():
If there is a passcode set, make a lockscreen appear first.
If userState is 0:
Display Normal in status box
If userState is -1:
Display Fatigued in status box
If userState is 1:
Display Stressed in status box

Graph Tab:

GraphTab:React.FC: ({userState, userSettings}) => JSX

Create states for graphData (any) and timeframeSelection (number)

getChartData: async () => void

Let timePeriod:number store the time period to look at in seconds
If timeframeSelection is 0:
timePeriod is set to 1 hour
Else If timeframeSelection is 1:
timePeriod is set to 1 day

```
Else if timeframeSelection is 2:
    timePeriod is set to 1 week
Else:
    Log an error to console
    End method
```

Let records:array be the output of fetchRecords(timePeriod)

```
Aggregate data:
    If timeframe is 1:
        Separate timestamps by nearest even hour
        Average the corresponding HRV values and timestamps and save in records
    If timeframe is 2:
        Separate timestamps into their corresponding day
        Average the corresponding HRV values and timestamps and save in records
```

Let labels:array be the local time representation of each timestamp. This is rounded to:
The nearest 5th minute when timeframe is 0,
The nearest Even hour when timeframe is 1,
The corresponding Day when timeframe is 2.

Let values:array be the corresponding HRV values of each record.

Let colors:array be the output of colorRecords(records).

```
Set chartData to a chartjs data object
labels should be labels
For the dataset:
    data should be values
    backgroundColor should be colors
    borderColor should be colors
```

```
colorRecords: async (records) => array
    Let colors:array be an empty array to store data point colors
```

```
Let baselineRecords:array be the output of fetchRecords(3 days (in seconds)).
Let baselineHRV:number be the mean of the HRV values stored in baselineRecords
```

```
For each record in records:
    If the record is fatigued or stressed (the record's HRV value is greater than 107 or greater than 115%
    or greater than userUpperThreshold of baselineHRV or less than 16 or less than 85% of baselineHRV or
    less than userLowerThreshold):
        Add the rgb value for a red color to the colors array
    Else if the record is greater than 108% of baselineHRV or less than 92% of baselineHRV:
        Add the rgb value for a yellow color to the colors array
    Else:
        Add the rgb value for a green color to the colors array
```

Return colors.

Return JSX object implementing the design of ()...

Advice Tab:
AdviceTab:React.FC () => JSX



Return JSX object implementing the design of () ...

Settings Tab:

SettingsTab:React.FC: ({userSettings}) => JSX

Return JSX object implementing the design of () ...:

Each text field edits its corresponding userSetting

Routing:

App:React.FC: () => JSX

Let userSettings: object be the a collection of useState<number> arrays for user upper and lower thresholds and passcode (userUpperThreshold, userLowerThreshold, passcode)

Let userState: object be the useState<number> array for the user state.

Let readjustError: number be the number of readjust errors from the last hour

hrvCallback: async (rawRecord: DataView) => void

storeRecord(rawRecord).
determineUserState().

determineUserState: async () => void

Let baselineRecords:array be the output of fetchRecords(3 days (in seconds)).

Let sampleRecords:array be the output of fetchRecords(3 minutes (in seconds)).

Let baselineHRV:number be the average of the baselineRecords' HRV values.

Let sampleHRV:number be the average of the sampleRecords' HRV values.

If the sample is stressed (the sampleHRV is greater than 107 or greater than 115% of baselineHRV or greater than userUpperThreshold):

Set userState to +1

Else if the sample is fatigued (the sampleHRV is or less than 16 or less than 85% of baselineHRV or less than userLowerThreshold):

Set userState to -1

Else:

Set userState to 0

storeRecord: async (rawRecord: DataView) => void

Parse rawRecord into timestamp:number and hrvMetric:number

Let currentDate: Date be the UTC date associated with timestamp

Open the corresponding HRV data file using currentDate (req...) for writing/appending

Let record:string be the formatted form of the record (req...)

If record is not already in the currentDate file:

Append the record to the file (end with newline).

errorCallback: async (rawError: DataView) => void

Parse rawError into errorCode:number

If errorCode is 1 (A readjust error):

Add one to readjustError

If readjustError is greater than 5:
Send a local notification to alert user to readjust their sensor

If readjustError is 1:
Set a timeout for an hour to set readjustError to 0.

Log errorCode to the console

Using useEffect (run once):
data_hook([hrvCallback, errorCallback])

Return JSX object with Routing to each page according to (), userState and userSettings should be passed as props to each component that requires it.

BLE Setup:

The device wearable acts as the BLE peripheral device

The device application acts as the BLE central device

HRV Service

UUID: x180F

HRV Characteristic

UUID: x2A19

Permissions: Read | Notify

Data Format: UUUUHHHH

U: Unix timestamp: uint32

H: HRV metric: float32

Error Characteristic

UUID: x2A1A

Permissions: Read | Notify

Data Format: Z

Z: Error Code: uint8

Request Characteristic

UUID: x2A1B

Permissions: Read | Write

Data Format: YYMD

Y: Year: uint16

M: Month: uint8

D: Day: uint8

TODO: make flowchart for between the devices, between each section of the app code, Arduino software priming directions, and app routing

Day 11: Design 4

Monday, July 31, 2023

6:24 PM

Goals:

- Consolidate all pseudocode into the Design Doc
 - Separate out function headers and variables
- Make software flowcharts
 - Between the devices
 - Between each section of the app code
 - Arduino software priming directions
 - app routing
- Time permitting, help with other sections
- Begin writing code for app, rtc init, finish DAQ code.